

MEMORANDUM FOR:

Ray

Tues or Thurs

*Would you call your
cleared? friend up at [redacted]
and advise him that two
people would like to see the
microdensitometric demonstration.*

*If he is not cleared please let
me know. These people are from*

[redacted]

14 Oct 64

(DATE)

[redacted]

FORM NO. 101 REPLACES FORM 10-101
1 AUG 54 WHICH MAY BE USED.

(47)

STATINTL

STATINTL

DECLASS REVIEW by NIMA/DOD

STATINTL

October 8, 1964

Attention: Raymond S.

STATINTL

You are cordially invited to inspect a [redacted]
System at our plant the week of 19 October 1964.

The complete system includes a high sensitivity microphotometer capable of reading optical density as high as 5.0 with a scanning area of 100 square microns, a linear measuring capability in two coordinates accurate to 1 part in 100,000, and digital equipment to record density as a function of X, Y coordinates on magnetic tape. Analog readout is also available. The instrument provides an automatic mode of operation with a programmed scan pattern. The instrument on display can also read density in three colors simultaneously.

Because you have previously expressed interest in [redacted] MicroSTATINTL densitometers, we have written you so that you may see such a system in actual operation.

If it is possible for you to visit, will you kindly telephone us to arrange a convenient day during the week. We shall be very pleased to show you this instrument as well as our new plant.

STATINTL

STATINTL

Very sincerely yours,

[redacted]
20 SEP 64
LOE:pc